

Intended use specified by the manufacturer:

Medical machine JETT PLASMA LIFT MEDICAL is a medical device that on the physical principle of succession spark discharges generated by a DC voltage is used for treatment of the skin, stops minor bleeding, removes small warts, haemangioma / benign so called benign mesenchymal carcinoma of blood vessels /and unwanted skin formations, venectasia, melasma, fibromas and keratosis.

The product is designed for surgical interventions above. It is a medical device for repeatable use marketed in non-sterile condition. The purpose of the device is the correction / removal of skin formations using spark discharges generated by a DC voltage.

Some of the above effects / clinical indications are not nature - it is a cosmetic defect. The following indications are specified as indications within the fields of cosmetics / not included in the conformity assessment process /:

- Skin Tags
- Small pedunculated growths
- Advanced veins
- Spider nevi
- Stretch marks
- Age spots

Indication:

- angioma senilis
- verrucae seborrhoicae
- verrucae plane
- angiokeratoma
- teleangiectasie
- lentigo
- fibroma molle
- keratoacanthoma
- conydomata acuminatum
- molluscum contagiosum
- verrucae vulgaris
- verrucae filiformes

- naveus capillaris
- naveus araneus
- keratosis actinica
- keratosis senilis
- keratosis seborrhoica
- melasma